

30th
Anniversary Celebration!

Your light at the end of tunnel



IMEXROJAS
IMPORTADORES - EXPORTADORES CIA. LTDA.

Ask your personal dealer rep for detailed pricing and product information. Sample is available upon request. Visit us at the following shows:

- AAE Show in San Francisco, April 6 - 8, 2016, Booth #807
- CDA Presents in Anaheim, May 12 -14, 2016, Booth #1284



NEW

Dia-PT™ File

Nickel Titanium Rotary File

FEATURES

Superior Cutting Efficiency

- Stronger and sharper blades can move deeper into the canal
- Triangular cross section reduces contact with the canal wall

High Flexibility

- Progressive taper design increases flexibility to reduce user fatigue
- High corrosion resistance with coated surface

Short Shaping Time

- Reduced shaping time versus traditional rotary technique
- Designed to instrument calcified or severely curved canals


Highly Integrable

- Designed to work interchangeably with your current technique and system
- Perfect matching GP and PP points: Dia-ProT™




Shaping

Finishing



{ Compatibility Chart }

- DX: Access File (19mm)
- D1: ProTaper® S1
- D2: ProTaper® S2
- D3: ProTaper® F1
- D4: ProTaper® F2
- D5: ProTaper® F3



DiaDent®

DiaDent Group International
www.diadent.com | 1.877.342.3368
Follow us on FACEBOOK

