

Elegan I Series CERAMIC BRACKET

Elegan I series brackets is the most successful classic aesthetic brackets in the markets. With the great aesthetic appearances and effective performance and the humbled price policy, Millions of pieces of brackets are sold monthly to every country of the world. Please don't miss this deal of best price/quality ratio!



Poly-crystalline aesthetics ceramic, CIM Technology.



IMEXROJAS

IMPORTADORES - EXPORTADORES CIA. LTDA.

Rank first in the world
Ceramic bracket series.



Advanced polishing technology, smooth surface.

Perfect radian design, and more comfortable for using.



Three groove base
Super good bonding.

Elegan I Series™ Roth System*

Patient Kits Description	Item Numbers	
	.018" Slot	.022" Slot
Upper/Lower 5-5 Hooks on Cuspid	E1-R18-01	E1-R22-01
Upper/Lower 5-5 Hooks on 3, 4, 5	E1-R18-02	E1-R22-02

▲ Maxillary						▼ Mandibulary							
Torq.	Ang.	R/L	.018" Slot	.022" Slot		Torq.	Ang.	R/L	.018" Slot	.022" Slot			
1 1	Anterior	+12	5	R	E1-R1811 E1-R1821	E1-R2211 E1-R2221	21 12	Lower Anterior	-1	0	U	E1-R1841	E1-R2241
2 2	2nd Insisor	+8	9	R	E1-R1812 E1-R1822	E1-R2212 E1-R2222	3 3	Cuspid with hook	-11	7	R	E1-R1843H E1-R1833H	E1-R2243H E1-R2233H
3 3	Cuspid with hook	-2	10	R	E1-R1813H E1-R1823H	E1-R2213H E1-R2223H	4 4	1" Bicuspid with hook	-17	0	R	E1-R1844H E1-R1834H	E1-R2244H E1-R2234H
54 45	Bicuspid with hook	-7	0	R	E1-R1814H E1-R1824H	E1-R2214H E1-R2224H	4 4	1" Bicuspid	-17	0	U	E1-R1844	E1-R2244
54 45	Bicuspid	-7	0	U	E1-R1814	E1-R2214	5 5	2" Bicuspid with hook	-22	0	R	E1-R1845H E1-R1835H	E1-R2245H E1-R2235H
							5 5	2" Bicuspid	-22	0	U	E1-R1845	E1-R2245

*The Creative Dental version of Roth is not claimed to be a duplication of any other, nor does Creative Dental imply that it is endorsed in any way by Dr. Roth.

Elegan I Series™ MBT System*

Patient Kits Description	Item Numbers	
	.018" Slot	.022" Slot
Upper/Lower 5-5 Hooks on Cuspid	E1-M18-01	E1-M22-01
Upper/Lower 5-5 Hooks on 3, 4, 5	E1-M18-02	E1-M22-02

▲ Maxillary						▼ Mandibulary							
Torq.	Ang.	R/L	.018" Slot	.022" Slot		Torq.	Ang.	R/L	.018" Slot	.022" Slot			
1 1	Anterior	+17	4	R	E1-M1811 E1-M1821	E1-M2211 E1-M2221	21 12	Lower Anterior	-6	0	U	E1-M1841	E1-M2241
2 2	2nd Insisor	+10	8	R	E1-M1812 E1-M1822	E1-M2212 E1-M2222	3 3	Cuspid with hook	-6	3	R	E1-M1843H E1-M1833H	E1-M2243H E1-M2233H
3 3	Cuspid with hook	-7	8	R	E1-M1813H E1-M1823H	E1-M2213H E1-M2223H	4 4	1" Bicuspid with hook	-12	2	R	E1-M1844H E1-M1834H	E1-M2244H E1-M2234H
54 45	Bicuspid with hook	-7	0	R	E1-M1814H E1-M1824H	E1-M2214H E1-M2224H	4 4	1" Bicuspid	-12	2	R	E1-M1844 E1-M1834	E1-M2244 E1-M2234
54 45	Bicuspid	-7	0	U	E1-M1814	E1-M2214	5 5	2" Bicuspid with hook	-17	2	R	E1-M1845H E1-M1835H	E1-M2245H E1-M2235H
							5 5	2" Bicuspid	-17	2	R	E1-M1845 E1-M1835	E1-M2245 E1-M2235

*The Creative Dental version of MBT is not claimed to be a duplication of any other, nor does Creative Dental imply that it is endorsed in any way by Drs. McLaughlin, Bennett, or Trevisi.