


Venus Series METAL BRACKET

- Laser welding technology.
- 17-4 stainless steel materials, ensures low friction and groove, never deformation.



The advanced technology of laser welding



3 times Polishing to ensure the bracket surface smooth and beautiful

Industrial best
0.002 precision error



80 mesh bottom
bracket design Mini



IMEXROJAS
IMPORTADORES - EXPORTADORES CIA. LTDA.

VENUS Series™ Roth System*

▲ Maxillary	Torq.	Ang.	R/L	.018" Slot	.022" Slot
1 1 Anterior	+12	5	R	VE-R1811 VE-R1821	VE-R2211 VE-R2221
2 2 2nd Insisor	+8	9	R	VE-R1812 VE-R1822	VE-R2212 VE-R2222
3 3 Cuspid with hook	-2	10	R	VE-R1813H VE-R1823H	VE-R2213H VE-R2223H
54 45 Bicuspid with hook	-7	0	R	VE-R1814H VE-R1824H	VE-R2214H VE-R2224H
54 45 Bicuspid	-7	0	U	VE-R1814	VE-R2214

*The Creative Dental version of Roth is not claimed to be a duplication of any other, nor does Creative Dental imply that it is endorsed in any way by Dr. Roth.

Patient Kits Description

Item Numbers

Upper/Lower 5-5 Hooks on Cuspids	VE-R18-01	.018" Slot	VE-R22-01	.022" Slot
Upper/Lower 5-5 Hooks on 3, 4, 5	VE-R18-02	.018" Slot	VE-R22-02	.022" Slot

VENUS Series™ MBT System*

▲ Maxillary	Torq.	Ang.	R/L	.018" Slot	.022" Slot
1 1 Anterior	+17	4	R	VE-M1811 VE-M1821	VE-M2211 VE-M2221
2 2 2nd Insisor	+10	8	R	VE-M1812 VE-M1822	VE-M2212 VE-M2222
3 3 Cuspid with hook	-7	8	R	VE-M1813H VE-M1823H	VE-M2213H VE-M2223H
54 45 Bicuspid with hook	-7	0	R	VE-M1814H VE-M1824H	VE-M2214H VE-M2224H
54 45 Bicuspid	-7	0	U	VE-M1814	VE-M2214

*The Creative Dental version of MBT is not claimed to be a duplication of any other, nor does Creative Dental imply that it is endorsed in any way by Drs. McLaughlin, Bennett, or Trevisi.

Patient Kits Description

Item Numbers

Upper/Lower 5-5 Hooks on Cuspids	VE-M18-01	.018" Slot	VE-M22-01	.022" Slot
Upper/Lower 5-5 Hooks on 3, 4, 5	VE-M18-02	.018" Slot	VE-M22-02	.022" Slot